

# Automated Sample Agitation Transportation System:

Cheap, Accessible Blood Test Automation

Dorjee Tenzing, Connor Platt, Gregory Ziegler, Mishari Alsabih, Alex Besanceney | Sponsored by Gamze Onuker  
ECE 129 | Stephen Petersen

ASATS


## Objective

Currently, an entire blood testing procedure takes two hours of repetitive, hands-on work, preventing nurses from providing direct care to patients. While this process has already been automated, the existing market solutions are expensive (~\$10k).

We have teamed up with PhD student Gamze Onuker to provide a cheaper, more customizable option which she will use in her innovative research to detect viral diseases earlier in their progression using magnetic bead-based purification.

## What our sponsor needs from us?

Customizability is key in Gamze's research. We have provided her sub-millimetre precision with speeds up to 300mm/s. We have also increased the number of test tube columns **1** she can use from 6 to 12, while simultaneously decreasing the number of user interactions from multiple to just one!




## Agitation Body Design

We added an additional axis of movement to agitate samples during blood testing. This design includes:

- A secondary body **2** that moves independently from a primary body **3**
- A slider to insert magnetic rods into plastic combs **4**

## A4988 Stepper Motor Driver

- Controls stepper motors with step and direction signals from the ESP32
- Each step turns the shaft 1.8°, but can be refined to 0.1125° with microstepping



## Power

- Buck regulator and SMPS regulate power from wall outlet
- Powers ESP32 microcontroller, stepper drivers and daughterboard





Figure: Power Supply and Motor Controller Motherboard

## How does the device work?

### Create a Protocol

User selects parameters for protocol



App


Scan Protocol  
onto  
Touchscreen



QR Scanner

### Select a Protocol

Daughterboard sends selected protocol data to motherboard via UART



Touchscreen

Begin Motor  
Operation




Scan me for demo!

### Executing a Protocol with Two Key Movements

Repeat the magnetize and agitation movements until the sample is purified.

**Magnetize:** bind sample to testing liquids



**Agitate:** remove contaminants by rapidly moving combs up and down

